
-ebro[®]

Elektronisches Kontaktthermometer
Electronic contact thermometer
Thermomètre de contact électronique
GFX 460 / GFX 460-G / GFX 460-B

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Table of Contents

Overview	5
The control elements	5
The display field	7
The measured value presentation	7
The heating display	7
Safety instructions	9
Unpacking/Delivery schedule	11
Delivery schedule	11
Installation	13
Device operation	15
Switching the device on/off	15
Changing the temperature desired value ..	15
Changing the measurement unit	17
Troubleshooting - what to do if	19
Cleaning and care	21
Calibration service	21
Disposal	21
Technical data	23
Declaration of Conformity	25

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
Overview

The control elements

The devices GFX 460, GFX 460 G and GFX 460-B are waterproof precision regulating thermometers for temperatures between -50°C to $+300^{\circ}\text{C}$ (-58°F to $+572^{\circ}\text{F}$). They are capable of controlling the heating of laboratory devices such as magnetic stirrers, heating stirrers etc.

Overshooting of the temperature is minimised by a fuzzy control. By means of this type of control the heating is switched off before reaching the desired value.

The GFX 460 device is fitted with a stainless steel probe (\varnothing 3 mm) and the GFX 460 G device with a glass probe (\varnothing 7 mm) for use with especially aggressive substances. Different Pt100 probes can be plugged on GFX 460-B.

Device elements:

- 1 Display field (LCD) with background lighting
- 2 Heating display (LED green/red)
- 3 Button ▼ (reduce desired temperature)
- 4 SET button
- 5 Temperature probe
- 6 connection to stirrer/current supply
- 7 Button ▲ ((increase desired temperature)

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

The display field

The display field (LCD) comprises the following display elements:

- 1 Tendency indicator ▲ / ▼
- 2 Current measured value
- 3 Desired value/text
- 4 Measurement unit (°C or °F)

The measured value presentation

During operation the device displays the temperature continuously in either °C or °F depending on the setting (see page 16).

The device updates the value displayed in one second cycles. When the temperature is rising or falling a tendency display appears to the left of the current temperature value.

The heating display

When the controls are set to heat, the device switches the heating contact by means of which a laboratory device is controlled to „on“. The colour of the light emitting diode under the display field now changes from green to red. Simultaneously »HEAT« appears in the display field below the current measured value.

The device switches the heating contact to „off“ again shortly before the measured value reaches the desired value. The light emitting diode changes back from red to green. The »HEAT« display fades out.

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Safety instructions

Do not use the device in potentially explosive areas. Mortal danger!

In order to avoid damage to the device, protect it from direct contact with water. Observe the protection class of the device (IP65)

Do not place the device in an environment that has temperatures exceeding 70°C. Potential danger of explosion of the rechargeable lithium battery!

Do not subject the device to microwave radiation. Potential danger of explosion of the rechargeable lithium battery!

Do not touch the probe with your bare hand immediately after measuring a high temperature. Risk of burns!

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Unpacking/Delivery schedule

Check that the contents are in perfect condition after receiving the device package. Also check that the contents of the device package comply with your order.

Delivery schedule

The delivery schedule comprises the following parts:

- GFX 460 Measuring Device with stainless steel probe, GFX 460-G with glass probe or GFX 460-B with plugable Pt100 probes.
- these operating instructions
- a calibration service card (see page 20)
- device fixture

In addition various accessories may be included in the delivery:

Description	Name
Mains adaptor	AX 410
Relay box	AX 400

Please contact us if you have reason to complain. You can find our contact address on the rear of the package.

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
Installation

The device needs a stable voltage supply of 9...12V. This is provided for many electrical laboratory devices by means of a 5 pin DIN plug socket.

For this reason the device has a connector cable with a DIN diode plug that also contains the heating control cable. This connector cable is designed for direct connection to the heating plates from the various leading manufacturers.

If needed you can find the plug configuration in the drawing on page 22.

You will need the relay box RB 1691 that is available as an accessory in order to connect heating plates or heating mantles that do not have a DIN plug socket.

- ◆ Positioning the device at the desired point.
- ◆ Put the DIN diode plug on the connector cable into the socket of the laboratory device to be used.

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


Device operation

Switching the device on/off

Current for the measuring device is supplied from the connected stirring device (optionally: by using the mains adapter provided). The measuring device starts up and is ready for use in seconds (self-testing) as soon as the stirring device is switched on or the mains adapter is plugged in.

The measuring device is switched off simultaneously as soon as the stirring device is switched off (or mains adapter is plugged out).

The display background lighting and the heating display (LED, see page 4) continue to light up for as long as the measuring device is supplied with current.

Changing the temperature desired value

- ◆ Press the »SET« button to change the temperature desired value.

The temperature desired value valid until now appears below the current measured value.

- ◆ Use the ▲ and ▼ buttons to change the temperature desired value.

One short press changes the desired value by 1°. For desired value changes of a greater range press and hold the button for a longer period.

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
Changing the measurement unit

You can set whether the device should measure the temperature in °Celsius or °Fahrenheit.

Call up the «User Menu» to make this setting. Now proceed as follows:

Prerequisite: the device is switched off (switch of stirring device or disconnect cable if necessary).

- ◆ Press and hold the »SET« button and then switch on current (switch on stirring device or connect cable or plug in optional mains adapter). Press and hold button until »USER« appears on display.

After a short wait »UNIT« appears on display.

- ◆ To change the measurement unit from °Celsius or °Fahrenheit or back again press »SET«.

The change is saved at once.

- ◆ To leave the «User Menu» switch off the device by switching off the connected stirring device, disconnecting the cable or pulling the mains adapter plug.

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Cleaning and maintenance

Clean the device with a slightly moist cloth.

Do not use any solvents such as acetone, because this can cause corrosion to the plastic.

Calibration service

The device has a high accuracy of measurement. In order to retain the accuracy, you should have the device calibrated annually. We offer a calibration service for this purpose.

◆ In order to be able to take advantage of this service, please fill in the enclosed service card and send it to ebro Electronic (address on the reverse side of these instructions), Fax +49 841 95478-80.

We will then collect your device and perform an annual calibration.

We will return the newly calibrated device to you within a week.

Disposal

If the device should become unfit for use, you must dispose of it in an appropriate and environmentally friendly manner.

Do not dispose of the device along with domestic waste. take it to an official collecting point for waste electrical and electronic equipment or return it to the manufacturer.

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
Technical data

Measurement	Temperature (°C/°F)
Probe GFX 460 /	GFX 460-G Pt1000
Probe GFX 460-B	Pt100 4-wire technologie
Measuring range	-50°C to +300°C -58°F to +572°F
Measuring precision	±1°C (±1°F)
Measurement interval	1second
Annul measuring interval	1°C, 1°F
Thermal time constant T99	
flowing water	
GFX 460 (Edelstahlmessfühler)	ca. 8s
GFX 460 G (Glasmessfühler)	ca. 8s
ruhendes Wasser	
GFX 460 (stainless steel probe)	approx. 10s
GFX 460 G (glass probe)	approx. 90s
Operating temperature	
	-5°C to +60°C (+23°F to +140°F)
Storage	
	-30°C to +70°C (-22°F to +158°F)
Control type	fuzzy control
System of protection	IP65
Current supply:	external 9-12V DC, 20mA
Safety features	Probe breakage recognition, probe terminal recognition, Safety circuit and control circuit to DIN 12878 classes 1 and 2
Other characteristics	Display background lighting LED indicator for control status
Connection	5 pin diode plug DIN 45322 see illustration left for pin assignment
Dimensions (L x W x H)	110 x 54 x 22 mm without probe
Housing material	Acid resistant polyamide
Weight	approx. 150 g

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ebro Electronic GmbH & Co. KG